

Veterinary ultrasound solution for Small and large animal pregnancy



- Designed for both small and large animal pregnancy — pigs, cows, sheep to dogs, cats and more
- High-quality images
- As easy to use, as it is easy to learn
- 340g (including the battery) light weight, small size and ergonomics designing
- More than 2.5 hours for continuous operation, ready anywhere, anytime with internal rechargeable battery and instant power-up
- Waterproof and dustproof probe, and can be easily cleaned by flowing water

Standard Configuration:

- One(l) main unit of M5 veterinary ultrasound system
- One(l) multi-frequency waterproof probe
- One(l) rechargeable Li Lon battery pack
- One(l) AC-adaptor (AC110V**240V(50/60Hz)
- One(l) wrist bandage
- One(l) carrying case
- Two(2) years factory warranty (one year for probe)

Optional:

- Car Lighter
- Video Goggles
- Additional rechargeable Li Lon battery pack



Specifications:

<ul style="list-style-type: none"> • Monitor: • Probe: • Frequency: • Depth: • Scan Angle: • Gray Shade: • Built-in Memory: • Common • OB 	<ul style="list-style-type: none"> • 3.5 inch high resolution color TFT LCD, 640 x 480 pixels • Multi-frequency mechanical sector • 2.5/3.5/5.0MHZ • Max 160mm • 78 degree • 256 gray levels • up to 128 frames • Distance, Area, Volume, Heart Rate • Pig-HL Sheep/Goat-USO • Bovine-BL, SL, HL Cat-HD, BD 	<ul style="list-style-type: none"> • AC-Adaptef: • Battery: • Weight: • Packing 	<ul style="list-style-type: none"> • Dog BD,GSO, CRI, HD, BD • 110V-240V, 50/60HZ • Internal Li lon battery pack • More than 2.5 hrsfor continuous • Fully charges batteries in 2hrs • 340g including the battery • 171mmx93mm5<62mm • Dimension 0.32x0.38x0.12m • N.W.(2.4KG), G.W.(3KG)
--	---	---	---